

## ENBIO S

VS

## ENBIO BEAUTY EDITION

**3 programs:** FAST in 7 min. (134°C),  
15 min. with drying (134°C),  
30 min. with drying (121°C)

**2 programs:** FAST in 7 min. (134°C),  
15 minutes with drying (134°C).

Ability to set up automatic tests:  
Helix/B&D and Vacuum.

\_\_\_\_\_

Ability to use "user identification"  
as a personalized signature

\_\_\_\_\_

Possibility of sterilization in the  
highest class - B (15 and 30 minutes)

Possibility of sterilization in the  
highest class - B (only in 15 min.)

**Design: medical.**  
Classic white cover, silver handle,  
modern interface.

**Design: fashion.**  
Gold elements and interface.  
Elegant color "secco" by Enbio.

enbio.