

BIOXTERIL[®]

Patented Solutions of Chlorine Dioxide (ClO₂)

- BIOXTERIL is suitable for **Osmo Steril** Trio y Pluri.
- BIOXTERIL P for **Osmo Steril** Mono.
- It provides constant disinfection and does not require continuous operator assistance.
- According to recent scientific studies m160s it is the most effective disinfectant and biocide on the market.



- It should be replaced every 30 days and only takes 2 minutes of time.
- It reduces the costs for water disinfection.

CODES

BIOXTERIL

Bicomponent system for 250ml of Chlorine Dioxide.

BIOXTERIL P

Bicomponent for 100ml of Chlorine Dioxide